

Electrical modules	Name	What	Input	Output	Supplier	Part no.	Location	Notes								
	T1	3 wire RTD			Omega	RTD-NPT-72-E	Top of top hat	TS-TE1 in the flow diagram								
	T2	4 wire RTD			Omega	PR-21D-3-100-A-1/8-1500-M12-1	Bottom of top hat	TS-TE2 in the flow diagram								
	T3	3 wire RTD			Omega	RTD-NPT-72-E	Top of water bath	TS-TE3 in the flow diagram								
	T4	3 wire RTD			Omega	RTD-NPT-72-E	Middle of water bath	TS-TE4 in the flow diagram								
	T5	3 wire RTD			Omega	RTD-NPT-72-E	Bottom of water bath	TS-TE5 in the flow diagram								
	PT1	pressure transducer	9-30 VDC (<10 mA)	0-5 VDC	Omega	PX319-300A5V	Hydraulic cart pressure									
	PT2	pressure transducer	9-30 VDC (<10 mA)	0-5 VDC	Omega	PX319-300A5V	High pressure accumulator									
	PT3	pressure transducer	9-30 VDC (<10 mA)	0-5 VDC	Omega	PX319-300A5V	Low pressure accumulator									
	PT4	pressure transducer	10-30 VDC	0-5 VDC	Setra	225 G250 PAD 42BB1	Chamber plumbing									
	PT5	pressure transducer	9-30 VDC (<10 mA)	0-5 VDC	Omega	PX319-300A5V	Top hat									
	PT6	dytran fast pressure transducer	0-5 VDC	0-5 VDC	Dytran	2005V	Top hat									
	Dytran supply	Dytran current supply		0-5 VDC, 1/4 A	Dytran	4112V										
	EV1	solenoid valve	120 VAC, 10 W		Omega	SV131	Compression valve									
	EV2	miniature solenoid valve	12, 24 VDC, 12 W max		Parker	009-0172-900	High pressure tweak on hydraulic system									
	EV3	miniature solenoid valve	12, 24 VDC, 12 W max		Parker	009-0172-900	Low pressure tweak on hydraulic system									
	Pump	miniature piston pump	120 VAC, 40 W		Clark	ET-150	Pump from low pressure accumulator to high pressure accumulator									
	H1	50 ohm resistive heater	50 W, 1 A		Cryocon	4039-012	Heater inside top hat of bubble chamber									
	Lakeshore controller	Temperature controller for heater	120 VAC	50 W, 1 A	Lakeshore	DRC-93C										
	Heat sense box	cRIO AI 6	1 A	1 A, 1 V	Custom		On top of Lakeshore									
	DI1	Digital input 1	TTL		Kings	KC-79-67 QD	Low voltage box									
	DI2	Digital input 2	TTL		Kings	KC-79-67 QD	Low voltage box									
	DI3	Digital input 3	TTL		Kings	KC-79-67 QD	Low voltage box									
	DI4	Digital input 4	TTL		Kings	KC-79-67 QD	Low voltage box									
	DO1	Digital output 1		TTL	Kings	KC-79-67 QD	Low voltage box									
	DO2	Digital output 2		TTL	Kings	KC-79-67 QD	Low voltage box									
	Ac1	Acoustic transducer 1		0-5 V	IUSB		Attached to bubble chamber									
	Mightex LEDs	Mightex LED driver Luxeon leds	120 VAC up to 1.4 A	24 V, 1 A DC, 3.0 A	Mightex Luxeon	SLC SA02 LXHL-LD3C	Inside Mightex box At chamber, next to bath	Typical 3 V drop across each								
	24 V supply 1		120 VAC, 60 W	24 V, 2.5 A	Rhino	PSC-24-060	Inside HV box									
	24 V supply 2		120 VAC, 60 W	24 V, 2.5 A	Rhino	PSC-24-060	Inside HV box									
	12 V supply		120 VAC, 15 W	12 V, 1.25 A	Rhino	PSC-12-015	Inside HV box									
	Relay 1	Solid state relay	120 VAC relay		Opto	120D3	Inside HV box									
	Relay 2	Solid state relay	120 VAC relay		Opto	120D3	Inside HV box									
	Relay 3	Solid state relay	24 V relay		Opto	DC60S3	Inside HV box									
	Relay 4	Solid state relay	24 V relay		Opto	DC60S3	Inside HV box									
	cRIO module	Sensor and digital i/o controller	24 V, 20 W MAX		NI	NI cRIO-9074	Inside LV box									
		DIO 5 V TTL cRIO module	TTL	TTL	NI	NI 9403	Connected to cRIO									
		AO +- 10 V		+-10 V	NI	NI 9263	Connected to cRIO									
		AI +- 20 mA	+- 20 mA		NI	NI 9203	Connected to cRIO									
		AI Universal	up to 60 V or 25 mA		NI	NI 9219	Connected to cRIO									
		AO 20 mA		20 mA	NI	NI 9265	Connected to cRIO									
		AI +-10 V	+- 10 V		NI	NI 9201	Connected to cRIO									
		AI RTD	RTD signals		NI	NI 9217	Connected to cRIO									
		DIO 5 V TTL cRIO module	TTL	TTL	NI	NI 9403	Connected to cRIO									
Manual valves	Name	Location	Supplier	Part number												
	MV1	Line from cart to bubble chamber	McMaster	45395K121												
	MV2	Needle valve behind EV2	McMaster	7836K26	Ideal Valve	5-2-11										
	MV3	Needle valve behind EV3	McMaster	7836K26												
	MV4	Fill line for high pressure accumulator	McMaster	45395K121												
	MV5	Pump port for high pressure accumulator	McMaster	45395K121												
	MV6	Fill line for low pressure accumulator	McMaster	45395K121												
	MV7	Pump port for low pressure accumulator	McMaster	45395K121												
	MV8	Fitting at top of high pressure accumulator	Parker	AD016B25T9A1												
	MV9	Fitting at top of low pressure accumulator	Parker	AD007B25T9A1												
	MV10	Needle valve on fluid side, before PT4	Swagelok	SS-4MG-VCR-MH												
	MV11	Fill line of inner volume	Swagelok	6LVV-DPFR4-P												
	MV12	Bleed out valve for inner volume	Swagelok	6LVV-DPFR4-P												
	MV13	Inlet line from hydraulic cart	McMaster	45395K121												
	MV14	Pump out port of hydraulic system on bubble chamber	McMaster	45395K112												
	MV15	Valve for neslab return	Swagelok	SS-4P4T4												
Miscellaneous	F1	Filter on fill line of inner volume	Pall	GLFPF3000VMM4												
	Hydraulic hose	Connection between hydraulic line and chamber	Parker	540N-6												
	</															

Source	Destination	Line	Cable wire color	Connector pins	Box wire color	Input pins	Notes
T1 - temperature sensor at top of top hat 3 wire RTD Omega RTD-NPT-72-E	Low voltage box, cRIO RTD module	+Excitation -Excitation Common	Red Black Black none	A B C D	Red Black Black none	RTD 1 RTD 2 RTD 3 none	TS-TE1 in flow diagram
T2 - temperature sensor at bottom of top hat 4 wire RTD Omega PR-21D-3-100-A-1/8-1500-M12-1	Lakeshore Temperature Controller Input A	Provided cable					TS-TE2 in flow diagram
Lakeshore Temperature Controller 7 pin output from back	Low voltage box, cRIO AI module	Voltage Out Voltage In 10 mV/K analog out (not active) Gnd for analog out Setpoint out Gnd for Voltage Out and Setpoint out Optional shield		A B C D E F H	Black White/Orange Red Black Green Brown White/Blue	AI 5   AI 9 AI 7 AI 9	
Lakeshore Temperature Controller Heater output	50 Ohm heater, via "Sense/Heat" box	Out Return	Coax lead Coax ground				Sense/Heat box contains a 1 ohm resistor to measure voltage across The heater has an isolated ground, the 1 ohm resistor comes after the heater. The box is at the Lakeshore ground
Sense/Heat box	cRIO, Heater Input, AI 6	V+ V-	Coax lead Coax ground		Red Black	AI 6 AI 9	
T3 - temperature sensor at top of bath 4 wire RTD Omega PR-21D-3-100-A-1/8-200-M12-1	Low voltage box, cRIO RTD module	+Excitation -Excitation +Signal -Signal	Brown White Blue Black	1-A 2-B 3-C 4-D	Red Black(Yellow tap) Black White	RTD 4 RTD 7 RTD 6 RTD 5	TS-TE3 in flow diagram
T4 - temperature sensor at mid bath 3 wire RTD Omega RTD-NPT-72-E	Low voltage box, cRIO RTD module	+Excitation -Excitation Common	Red Black Black none	A B C D	Red Black Black none	RTD 9 RTD 10 RTD 11 none	TS-TE4 in flow diagram
T5 - temperature sensor at bottom of bath 3 wire RTD Omega PR-21D-3-100-A-1/8-200-M12-1	Low voltage box, cRIO RTD module	+Excitation -Excitation +Signal -Signal	Brown White Blue Black	1-A 2-B 3-C 4-D	Red Black(Yellow tap) Black White	RTD 4 RTD 7 RTD 6 RTD 5	TS-TE5 in flow diagram
PT1 - pressure at hydraulic system Omega PX319-300A5V	Low voltage box, cRIO AI module	+Excitation -Excitation Sgn	Red Green Black none	A B C D	Red Green Black none	+12 V junction rail -12 V junction rail AI 0 none	-12 V junction rail is tied to AI 9 (common ground on AI module)
PT2 - high pressure accumulator Omega PX319-300A5V	Low voltage box, cRIO AI module	+Excitation -Excitation Sgn	Red Green Black none	A B C D	Red Green Black none	+12 V junction rail -12 V junction rail AI 1 none	-12 V junction rail is tied to AI 9 (common ground on AI module)
PT3 - low pressure accumulator Omega PX319-300A5V	Low voltage box, cRIO AI module	+Excitation -Excitation Sgn	Red Green Black none	A B C D	Red Green Black none	+12 V junction rail -12 V junction rail AI 2 none	-12 V junction rail is tied to AI 9 (common ground on AI module)
PT4 - bubble chamber pressure Setra 225 G250 PAD 42BB1	Low voltage box, cRIO AI module	+Excitation +Output -Output -Excitation	Red Black Green White	A B C D	Red White Green Black	+12 V junction rail AI 3 AI 9 -12 V junction rail	-12 V junction rail is tied to AI 9 (common ground on AI module)
PT5 - top hat pressure Omega PX319-300A5V	Low voltage box, cRIO AI module	+Excitation -Excitation Sgn	Red White Green none	A B C D	Red Green Black none	+12 V junction rail -12 V junction rail AI 4 none	-12 V junction rail is tied to AI 9 (common ground on AI module)
PT6 - Dytran fast pressure transducer Dytran 2005V	Dytran Current source	Provided cable					
Mightex box - trigger	Low voltage box, cRIO DIO module	+Cam 0 -Cam 0 +Cam 1 -Cam 1	Red Orange Black Brown	A - A B - B H - C J - D	Red Orange Green Black	DI 0 DI 9 (cmn gnd) DI 1 DI 10 (cmn gnd)	
Mightex box - LED driver	LED array	+LED 0 - LED0 +LED 1 -LED 1	Brown Red Black Orange	A B H J		LED 0 + LED 0 - LED 1+ LED 1-	
LED array Cam0 Cam1	From Mightex box From Mightex via RJ45 cable From Mightex via RJ45 cable						
Acoustic transducer	Acoustic breakout box	Vcc Vee Vgnd Sgn+ Sgn-	Red Green Black Blue White	A C B F E	Green Red Black Black White	Internal 3 Internal 4 Internal 5 Internal 1 Internal 2	Red White Black Green Orange
Trigger signals from PC and Labview	DI 1-4, to cRIO DIO module BNC inputs	+TTL Gnd			Red Black	DI 16-19 DI cmn gnd	
	DO 1-4, to cRIO DIO module BNC outputs	+TTL Gnd			Red Black	DO 24-27 DO cmn gnd	
EV1 - compression valve	High voltage box	+120 AC Neutral	Black Black	Conduit		+120 AC junction rail Gnd junction rail	Junction rail powered by relay 2, which is activated by DO 0 of cRIO
EV2 - high pressure tweak Parker 009-0172-900	High voltage box	+24 V -24 V	Yellow Yellow	Conduit		+24 V junction rail -24 V junction rail	Junction rail powered by relay 3, which is activated by DO 1 of cRIO
EV3 - low pressure tweak Parker 009-0172-900	High voltage box	+24 V -24 V	Yellow Yellow	Conduit		+24 V junction rail -24 V junction rail	Junction rail powered by relay 4, which is activated by DO 2 of cRIO
Pump Clark ET-150 pump	High voltage box	+120 AC Neutral Gnd	Black Gray Green	Conduit		120 VAC junction rail Neutral junction rail Gnd junction rail	Junction rail powered by relay 1, which is activated by DO 3 of cRIO
High voltage box Power and relay lines to control system	Low voltage box Voltage to cRIO, pressure transducers Relay lines to cRIO DIO module	+24 V -24 V +12 V -12 V Relay 1 Relay 2 Relay 3 Relay 4 Relay Gnd	Red (Yellow tape) Black (Yellow tape) Red (Blue tape) Black (Blue tape) White/Brown White/Yellow White/Orange White/Blue White/Gray	Conduit		cRIO +24 input cRIO -24 input +12 V junction rail -12 V junction rail DO 3 (Pin 4) DO 0 (Pin 1) DO 1 (Pin 2) DO 2 (Pin 3) DO 9 (Pin 10)	power for pressure transducers
AC Power supply attached to high voltage box	High voltage box	+120 AC Neutral Gnd	Black Gray Green	Conduit		120 VAC junction rail Neutral junction rail Gnd junction rail	Provides power for most of system

	Old T3							
	Low voltage box	+Excitation	Red	A	Red	RTD 5		
		-Excitation	Black	B	Black	RTD 6		
		Common	Black	C	Black	RTD 7		
			none	D	none	none		
	Old T5							
	T5 - temperature	Low voltage box	+Excitation	Red	A	Red	RTD 13	TS-TE5 in flow c
	3 wire RTD		-Excitation	Black	B	Black	RTD 14	
	Omega PR-21D-3-100-A-1/8-20	Common	Black	C	Black	RTD 15		
				none	D	none	none	